



**Health Fast Surgical**  
Surgical Manufacturer Company


# Orthodontics




## **DENTAL INSTRUMENTS CATALOGUE**




Drahtbiegezaugen  
Wire Bending Pliers  
Pinze per piegare fili  
Alicates para doblar alambres  
Pincers à pliers les fils



DHDP-118  
130 mm




DI-DP-119  
130 mm




DI-DP-121  
130 mm


**Tungsten-Carbide Inserts**  
Mit Hartmetalleinsätzen  
Plaquettes en carbure de tungstène  
Plaquetas de carburo de tungsteno  
Inserti in carburo di tungsteno



DI-LX-84  
125 mm



DI-LX-85  
125 mm




DI-LX-87  
125 mm




Drahtbiegezaugen  
Wire Bending Pliers  
Pinze per piegare fili  
Alicates para doblar alambres  
Pincas à pliers les fils

TC




DI-DP-114  
130 mm



Drahtbiegezaugen  
Wire Bending Pliers  
Pinze per piegare fili  
Alicates para doblar alambres  
Pincas à pliers les fils

Mod. Medium  
Mod. Medio




ADERER  
DI-DP-116  
125 mm




TC

**Tungsten-Carbide Inserts**  
Mit Hartmetalleinsätzen  
Plaquettes en carbure de tungstène  
Plaquetas de carburo de tungsteno  
Inserti in carburo di tungsteno




DI-LX-101  
185 mm

TC




WEINGART  
DI-DP-135  
140 mm



WEINGART  
DI-DP-134  
140 mm



Drahtschneidezangen  
Wire cutting pliers  
Tronchesi per fili metallici  
Alicates para cortar alambre  
Pinces à couper les fils TC



DI-PH-58  
150 mm




SMITH  
DI-FG-161  
160 mm




SMITH  
DI-FG-160  
160 mm

**Tungsten-Carbide Inserts**  
Mit Hartmetalleinsätzen  
Plaquettes en carbure de tungstène  
Plaquetas de carburo de tungsteno  
Inserti in carburo di tungsteno



DI-LX-102  
145 mm




○ Max Ø 2.0 mm (.078")


Für weichen Draht  
For soft Wire  
Per filo leggero  
Para alambre blando  
Pour fil léger

○ Max Ø 1.5 mm (.059")

Für harten Draht  
For hard Wire  
Per filo duro  
Para alambre duro  
Pour fil dur



DI-LX-103  
170 mm



○ Max Ø 2.0 mm (.078")

Für weichen Draht  
For soft Wire  
Per filo leggero  
Para alambre blando  
Pour fil léger

○ Max Ø 1.5 mm (.059")



Für harten Draht  
For hard Wire  
Per filo duro  
Para alambre duro  
Pour fil dur







Scheren  
Scissors  
Forbici  
Tijeras  
Ciseaux

Scheren für dünne Stahl-und Edelmetallbleche - Scissors for thin steel and precious metal sheets  
Forbici per lamine sottili di acciaio e metalli preziosi - Tijeras para láminas delgadas de acero y metales preciosos  
Ciseaux pour couper les feuilles de métaux



Spitz  
Sharp  
Acute  
Aguda  
Pointu



BEEBE  
DI-DP-91/ DI-DP-92  
110 mm



Stumpf  
Blunt  
Smusse  
Roma  
Mousse


BEEBE  
DI-DP-102/ DI-DP-103  
110 mm

Scheren für dünne Stahl-und Edelmetallbleche - Scissors for thin steel and precious metal sheets  
Forbici per lamine sottili di acciaio e metalli preziosi - Tijeras para láminas delgadas de acero y metales preciosos  
Ciseaux pour couper les feuilles de métaux



Draht-und Blechschere gezahnt  
Wire cutting scissor serrated  
Forbice universale dentellata  
Tijeras para cortar hilo  
de acero dentadas  
Coupe-fils dentelé

DI-DP-90  
120 mm



**Tungsten-Carbide Inserts**  
Mit Hartmetalleinsätzen  
Plaquettes en carbure de tungstène  
Plaquetas de carburo de tungsteno  
Inserti in carburo di tungsteno



o Max Ø 1.2 mm (.048")

Für weichenDraht  
For soft Wire  
Per filo leggero  
Para alambre blando  
Pour fil léger




DI-DP-120  
120 mm



Separierzangen  
Separating Pliers  
Pinze tendielastici  
Pinza estira-elásticos  
Pince à separer



Separierzange, zum Einbringen von Gummi-Separierligaturen  
Separating forceps for inserting elastic separators  
Pinze per legature per separatori elastici  
Pinza para ligaduras estira-elásticos  
Pince à separer les ligatures




Ligaturenspannzange  
Ligature Tying Plier  
Pinza per legature  
Pinza para ligaduras  
Pince à serrer les ligatures



Bracketset Pinzette  
Bonding Bracket Forceps  
Pinzette per brackets  
Pinza fina para ligar brackets  
Fine pince à lier des boîtiers

5.0  
4.0  
0



**Aluminium**  
Aluminium  
Aluminium  
Aluminio  
Alluminio

DI-DP-100  
170 mm

Bracketpositioner  
Bracket Positioning Gauge  
Posizionatore per attacchi  
Posicionador de Brackets  
Positionateur des Brackets



Gezant  
Serrated  
Dentellato  
Dentado  
Dentelé

MERSHON  
DI-DP-50  
140 mm

Verbandschieber  
Band Pusher  
Spingibande  
Empujador de bandas  
Enforcee-bague

Utility Instrumente  
Utility Instruments  
Strumenti di utilità  
Instrumentos para utilidad  
Instruments pour utility


Ligaturen Instrumente  
Ligating Instruments  
Strumenti per legature  
Instrumentos para ligaduras  
Instruments pour ligatures

DI-DP-100  
170 mm





SHURE  
DI-DP-124  
175 mm


2,5 x 1,0 mm




DI-DP-54  
175 mm



DI-DP-125  
155 mm



DI-PF-1  
170 mm




DI-PF-2  
170 mm




DI-PF-3  
160 mm




Ligature Instruments  
Ligating Instruments  
Strumenti per legature  
Instrumentos para ligaduras  
Instruments pour ligatures



PEAN-MURPHY  
DI-BF-22  
130 mm



HALSTED-MOSQUITO  
DI-P-96/ DI-P-97  
120 mm



DI-COLLEGE  
150 mm



MILLER  
DI-DR-78  
145 mm

Bracketset Pinzette - Forceps for positioning brackets  
Pinzette per riposizionare gli attacchi diretti  
(posizionamento dei bracket)  
Pinza para colocar brackets - Précelle pour positionner les brackets

Bracket Gauge  
Stainless Steel



DI-DP-132  
130 mm




DI-DP-133  
130 mm



DI-DP-52  
130 mm




DI-DP-126  
110 mm



DI-DP-127  
125 mm



DI-DP-128  
145 mm




DI-DP-129  
130 mm




Messinstrumente  
Measuring Instruments  
Strumenti per misura  
Instrumentos de medición  
Instruments à mesurer

Dental-Schieblehre für Praxis und Labor  
Gauge Caliper for practice and laboratory  
Calibri dentali per studio e laboratorio  
Calibrador dental para práctica y el laboratorio  
Gauge calibreur pour practice et laboratoire




0 - 10 mm  
Messbereich  
Measuring range  
Campo di misura  
Rango de medición  
Champ de mesure

IWANSON  
DI-DP-721  
100 mm




0 - 10 mm  
Messbereich  
Measuring range  
Campo di misura  
Rango de medición  
Champ de mesure


IWANSON  
DI-DP-720  
100 mm




DI-DP-723  
80 mm



DI-DP-724  
80 mm



ZIELINSKY  
DI-DP-725  
110 mm




BOLEY  
DI-DP-730  
140 mm

0 - 100 mm  
Messbereich  
Measuring range  
Campo di misura  
Rango de medición  
Champ de mesure






Zubehör  
Accessories  
Accessori  
Accesorios  
Accessoires



Zangenständer  
Plier Stand  
Supporto per pinze  
Soportes para alicates  
Supports de pinces

**Stainless Steel**  
Rostfreier Stahl  
Acier inoxydable  
Acero inoxidable  
Acciaio inossidabile



DI-DP-140